

For Research Use Only

Claudin 18 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60782-2-PBS



Basic Information

Catalog Number:

60782-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34049

GenBank Accession Number:

NM_001002026

GeneID (NCBI):

51208

UNIPROT ID:

P56856-2

Full Name:

claudin 18

Calculated MW:

28 kDa

Purification Method:

Protein G purification

CloneNo.:

1A12E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51118-1, Claudin 18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60782-1-PBS. Detection antibody: 60782-2-PBS. Standard:Ag34049. Range: 6.25-100 ng/mL.



Cytometric bead array standard curve of MP51118-2, Claudin 18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60782-3-PBS. Detection antibody: 60782-2-PBS. Standard:Ag34049. Range: 3.125-100 ng/mL.