

For Research Use Only

NHSL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60784-4-PBS



Basic Information

Catalog Number:

60784-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG19497

GenBank Accession Number:

BC136756

GeneID (NCBI):

340527

UNIPROT ID:

Q5HYW2

Full Name:

NHS-like 2

Calculated MW:

709 aa, 76 kDa

Purification Method:

Protein A purification

CloneNo.:

1H10E1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51122-3, NHSL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60784-1-PBS. Detection antibody: 60784-4-PBS. Standard:Ag19497. Range: 6.25-100 ng/mL