

For Research Use Only

# LAMA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60790-1-PBS



## Basic Information

Catalog Number:	60790-1-PBS	GenBank Accession Number:	NM_005559	Purification Method:	Protein G purification
Size:	1 mg/ml	GeneID (NCBI):	284217	CloneNo.:	3E6D2
Source:	Mouse	UNIPROT ID:	P25391		
Isotype:	IgG1	Full Name:	laminin, alpha 1		
Immunogen Catalog Number:	AG29326	Calculated MW:	337 kDa		

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51131-1, LAMA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60790-1-PBS. Detection antibody: 60790-2-PBS. Standard:Ag29326. Range: 0.098-100 ng/mL.