

For Research Use Only

TEK Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60793-2-PBS



Basic Information

Catalog Number:

60793-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13523

GenBank Accession Number:

BC035514

GeneID (NCBI):

7010

UNIPROT ID:

Q02763

Full Name:

TEK tyrosine kinase, endothelial

Calculated MW:

1124 aa, 126 kDa

Purification Method:

Protein G purification

CloneNo.:

2A1C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51137-1, TEK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60793-1-PBS. Detection antibody: 60793-2-PBS. Standard:Ag13523. Range: 0.781-200 ng/mL.