

For Research Use Only

MERTK Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60794-1-PBS



Basic Information

Catalog Number:

60794-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG17401

GenBank Accession Number:

BC114918

GeneID (NCBI):

10461

UNIPROT ID:

Q12866

Full Name:

c-met proto-oncogene tyrosine kinase

Calculated MW:

999 aa, 110 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A6D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51138-1, MERTK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60794-1-PBS. Detection antibody: 60794-2-PBS. Standard:Ag17401. Range: 0.098-25 ng/mL