

For Research Use Only

EPB41 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60801-1-PBS



Basic Information

Catalog Number:

60801-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4365

GenBank Accession Number:

BC039079

GeneID (NCBI):

2035

UNIPROT ID:

P11171

Full Name:

erythrocyte membrane protein band 4.1 (elliptocytosis 1, RH-linked)

Calculated MW:

720 aa, 81 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C7C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51150-1, EPB41 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60801-1-PBS. Detection antibody: 60801-2-PBS. Standard:Ag4365. Range: 3.125-100 ng/mL.



Biolayer interferometry (BLI) kinetic assays of 60801-1-Ig against Human EPB41 were performed. The affinity constant is 0.609 nM.