

For Research Use Only

ADRA2B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60802-2-PBS



Basic Information

Catalog Number: 60802-2-PBS	GenBank Accession Number: NM_000682.7	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 151	CloneNo.: 2H4H1
Source: Mouse	UNIPROT ID: P18089	
Isotype: IgG1	Full Name: adrenergic, alpha-2B-, receptor	
Immunogen Catalog Number: AG37590	Calculated MW: 50kDa, 450aa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51152-1, ADRA2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60802-1-PBS. Detection antibody: 60802-2-PBS. Standard:Ag37590. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP51152-2, ADRA2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60802-3-PBS. Detection antibody: 60802-2-PBS. Standard:Ag37590. Range: 1.563-50 ng/mL.