

For Research Use Only

ADRA2B Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60802-3-PBS



Basic Information

Catalog Number:

60802-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG37590

GenBank Accession Number:

NM_000682.7

GeneID (NCBI):

151

UNIPROT ID:

P18089

Full Name:

adrenergic, alpha-2B-, receptor

Calculated MW:

50kDa, 450aa

Purification Method:

Protein A Magarose purification

CloneNo.:

1H10F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51152-2, ADRA2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60802-3-PBS. Detection antibody: 60802-2-PBS. Standard:Ag37590. Range: 1.563-50 ng/mL



Cytometric bead array standard curve of MP51152-4, ADRA2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60802-1-PBS. Detection antibody: 60802-3-PBS. Standard:Ag37590. Range: 3.125-100 ng/mL