

For Research Use Only

LEF1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60803-1-PBS



Basic Information

Catalog Number: 60803-1-PBS	GenBank Accession Number: BC050632	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 51176	CloneNo.: 1B4F1
Source: Mouse	ENSEMBL Gene ID: ENSG00000138795	
Isotype: IgG1	UNIPROT ID: Q9UJU2	
Immunogen Catalog Number: AG29629	Full Name: Lymphoid enhancer-binding factor 1	
	Calculated MW: 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51153-1, LEF1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60803-1-PBS. Detection antibody: 60803-2-PBS. Standard:Ag29629. Range: 1.563-25 ng/mL