

For Research Use Only

ABCA5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60804-2-PBS



Basic Information

Catalog Number:

60804-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG21956

GenBank Accession Number:

BC054480

GeneID (NCBI):

23461

UNIPROT ID:

Q8WWZ7

Full Name:

ATP-binding cassette, sub-family A

(ABC1), member 5

Calculated MW:

1642 aa, 187 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A1B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51155-1, ABCA5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60804-2-PBS. Detection antibody: 60804-3-PBS. Standard:Ag21956. Range: 0.391-12.5 ng/mL