

For Research Use Only

ABCA5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60804-3-PBS



Basic Information

Catalog Number: 60804-3-PBS	GenBank Accession Number: BC054480	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 23461	CloneNo.: 1B12F3
Source: Mouse	UNIPROT ID: Q8WWZ7	
Isotype: IgG1	Full Name: ATP-binding cassette, sub-family A (ABC1), member 5	
Immunogen Catalog Number: AG21956	Calculated MW: 1642 aa, 187 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51155-1, ABCA5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60804-2-PBS. Detection antibody: 60804-3-PBS. Standard:Ag21956. Range: 0.391-12.5 ng/mL