

For Research Use Only

# AOC2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60806-2-PBS



## Basic Information

Catalog Number:

60806-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26965

GenBank Accession Number:

BC142641

GeneID (NCBI):

314

UNIPROT ID:

O75106

Full Name:

amine oxidase, copper containing 2 (retina-specific)

Calculated MW:

756 aa, 84 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

6F2B7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51161-1, AOC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60806-1-PBS. Detection antibody: 60806-2-PBS. Standard:Ag26965. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP51161-2, AOC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60806-3-PBS. Detection antibody: 60806-2-PBS. Standard:Ag26965. Range: 6.25-200 ng/mL