

For Research Use Only

PHF11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60808-1-PBS



Basic Information

Catalog Number:

60808-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG18041

GenBank Accession Number:

BC017212

GeneID (NCBI):

51131

UNIPROT ID:

Q9UIL8

Full Name:

PHD finger protein 11

Calculated MW:

34 kDa

Purification Method:

Protein A purification

CloneNo.:

1D4D5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51172-1, PHF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60808-1-PBS. Detection antibody: 60808-2-PBS. Standard:Ag18041. Range: 0.195-200 ng/mL



Cytometric bead array standard curve of MP51172-2, PHF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60808-1-PBS. Detection antibody: 60808-3-PBS. Standard:Ag18041. Range: 1.563-200 ng/mL