

For Research Use Only

CUZD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60809-2-PBS



Basic Information

Catalog Number:

60809-2-PBS

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG24335

GenBank Accession Number:

BC136755

GeneID (NCBI):

50624

UNIPROT ID:

Q86UP6

Full Name:

CUB and zona pellucida-like domains

1

Calculated MW:

68 kDa

Purification Method:

Caprylic acid/ammonium sulfate precipitation

CloneNo.:

1A5C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51173-1, CUZD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60809-1-PBS. Detection antibody: 60809-2-PBS.



Cytometric bead array standard curve of MP51173-1, CUZD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60809-1-PBS. Detection antibody: 60809-2-PBS. Standard: Ag24335. Range: 0.781-100 ng/mL.