

For Research Use Only

# CUZD1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60809-3-PBS



## Basic Information

Catalog Number:  
60809-3-PBS

Source:  
Mouse

Isotype:  
IgM

Immunogen Catalog Number:  
AG24335

GenBank Accession Number:  
BC136755

GeneID (NCBI):  
50624

UNIPROT ID:  
Q86UP6

Full Name:  
CUB and zona pellucida-like domains  
1

Calculated MW:  
68 kDa

Purification Method:  
Caprylic acid/ammonium sulfate  
precipitation

CloneNo.:  
1D9D5

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51173-2, CUZD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60809-1-PBS. Detection antibody: 60809-3-PBS.



Cytometric bead array standard curve of MP51173-2, CUZD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60809-1-PBS. Detection antibody: 60809-3-PBS. Standard: Ag24335. Range: 0.781-100 ng/mL.