

For Research Use Only

TMEM168 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60811-2-PBS



Basic Information

Catalog Number:

60811-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18337

GenBank Accession Number:

BC030804

GeneID (NCBI):

64418

UNIPROT ID:

Q9H0V1

Full Name:

transmembrane protein 168

Calculated MW:

697 aa, 80 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C12F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51180-1, TMEM168 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60811-1-PBS. Detection antibody: 60811-2-PBS. Standard:Ag18337. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51180-2, TMEM168 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60811-3-PBS. Detection antibody: 60811-2-PBS. Standard:Ag18337. Range: 0.391-100 ng/mL