

For Research Use Only

CCR3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60815-2-PBS



Basic Information

Catalog Number:

60815-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27235

GenBank Accession Number:

BC110297

GeneID (NCBI):

1232

UNIPROT ID:

P51677

Full Name:

chemokine (C-C motif) receptor 3

Calculated MW:

355 aa, 41 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1H9B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51186-1, CCR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60815-1-PBS. Detection antibody: 60815-2-PBS. Standard:Ag27235. Range: 6.25-200 ng/mL



Cytometric bead array standard curve of MP51186-4, CCR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60815-5-PBS. Detection antibody: 60815-2-PBS. Standard:Ag27235. Range: 0.195-12.5 ng/mL