

For Research Use Only

# CCR3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60815-4-PBS



## Basic Information

Catalog Number: 60815-4-PBS	GenBank Accession Number: BC110297	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 1232	CloneNo.: 2B6A4
Source: Mouse	UNIPROT ID: P51677	
Isotype: IgG1	Full Name: chemokine (C-C motif) receptor 3	
Immunogen Catalog Number: AG27235	Calculated MW: 355 aa, 41 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51186-3, CCR3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60815-4-PBS. Detection antibody: 60815-5-PBS. Standard:Ag27235. Range: 3.125-200 ng/mL.