

For Research Use Only

# RAB26 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60816-3-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60816-3-PBS	BC066913	Protein G purification
Size:	GeneID (NCBI):	CloneNo.:
1 mg/ml	25837	2B3D7
Source:	UNIPROT ID:	
Mouse	Q9ULW5	
Isotype:	Full Name:	
IgG1	RAB26, member RAS oncogene family	
Immunogen Catalog Number:	Calculated MW:	
AG6463	28 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51188-2, RAB26 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60816-1-PBS. Detection antibody: 60816-3-PBS. Standard:Ag6463. Range: 0.781-50 ng/mL