

For Research Use Only

# GYPC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60820-1-PBS



## Basic Information

Catalog Number:

60820-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG12825

GenBank Accession Number:

BC106051

GeneID (NCBI):

2995

UNIPROT ID:

P04921

Full Name:

glycophorin C (Gerbich blood group)

Calculated MW:

128 aa, 14 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E12A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51194-1, GYPC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60820-1-PBS. Detection antibody: 60820-2-PBS. Standard:Ag12825. Range: 0.781-100 ng/mL.