

For Research Use Only

# CXCL7/PPBP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60823-1-PBS



## Basic Information

Catalog Number:

60823-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4147

GenBank Accession Number:

BC028217

GeneID (NCBI):

5473

UNIPROT ID:

P02775

Full Name:

pro-platelet basic protein (chemokine (C-X-C motif) ligand 7)

Calculated MW:

128 aa, 14 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3A2B10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51217-1, PPBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60823-1-PBS. Detection antibody: 60823-2-PBS. Standard:Ag4147. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51217-2, PPBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60823-1-PBS. Detection antibody: 60823-3-PBS. Standard:Ag4147. Range: 0.391-100 ng/mL