

For Research Use Only

TMEM176B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60825-1-PBS



Basic Information

Catalog Number:

60825-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13918

GenBank Accession Number:

BC004358

GeneID (NCBI):

28959

UNIPROT ID:

Q3YBM2

Full Name:

transmembrane protein 176B

Calculated MW:

270 aa, 29 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D1A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51220-1, TMEM176B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60825-1-PBS. Detection antibody: 60825-2-PBS. Standard:Ag13918. Range: 0.098-100 ng/mL