

For Research Use Only

VWA1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60826-4-PBS



Basic Information

Catalog Number:

60826-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG37064

GenBank Accession Number:

BC059409

GeneID (NCBI):

64856

UNIPROT ID:

Q6PCB0

Full Name:

von Willebrand factor A domain containing 1

Calculated MW:

47 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F1C9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51222-3, VWA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60826-1-PBS. Detection antibody: 60826-4-PBS. Standard:Ag37064. Range: 0.195-100 ng/mL.