

For Research Use Only

ATP8A1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60828-2-PBS



Basic Information

Catalog Number:

60828-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16071

GenBank Accession Number:

BC109317

GeneID (NCBI):

10396

UNIPROT ID:

Q9Y2Q0

Full Name:

ATPase, aminophospholipid transporter (APLT), class I, type 8A, member 1

Calculated MW:

1164 aa, 131 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B2C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51223-1, ATP8A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60828-1-PBS. Detection antibody: 60828-2-PBS. Standard:Ag16071. Range: 0.781-100 ng/mL.



Biolayer interferometry (BLI) kinetic assay of 60828-2-Ig against Human ATP8A1 was performed. The affinity constant is 0.447 nM...