

For Research Use Only

NHE8 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60829-1-PBS



Basic Information

Catalog Number:

60829-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG13102

GenBank Accession Number:

BC112213

GeneID (NCBI):

23315

UNIPROT ID:

Q9Y2E8

Full Name:

solute carrier family 9
(sodium/hydrogen exchanger),
member 8

Calculated MW:

581 aa, 65 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B3E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51224-1, NHE8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60829-1-PBS. Detection antibody: 60829-2-PBS. Standard:Ag13102. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP51224-2, NHE8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60829-1-PBS. Detection antibody: 60829-3-PBS. Standard:Ag13102. Range: 0.098-100 ng/mL