

For Research Use Only

RPGR Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60830-3-PBS



Basic Information

Catalog Number: 60830-3-PBS	GenBank Accession Number: BC031624	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 6103	CloneNo.: 1G2E3
Source: Mouse	UNIPROT ID: Q92834	
Isotype: IgG1	Full Name: retinitis pigmentosa GTPase regulator	
Immunogen Catalog Number: AG27524	Calculated MW: 113 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51225-2, RPGR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60830-1-PBS. Detection antibody: 60830-3-PBS. Standard:Ag27524. Range: 0.195-100 ng/mL.