

For Research Use Only

SFTPD Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60831-1-PBS



Basic Information

Catalog Number:

60831-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2461

GenBank Accession Number:

BC022318

GeneID (NCBI):

6441

UNIPROT ID:

P35247

Full Name:

surfactant protein D

Calculated MW:

375 aa, 38 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F11C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51227-1, SFTPD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60831-1-PBS. Detection antibody: 60831-2-PBS. Standard:Ag2461. Range: 0.098-50 ng/mL