

For Research Use Only

DEPDC5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60832-2-PBS



Basic Information

Catalog Number:

60832-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG30176

GenBank Accession Number:

BC146766

GeneID (NCBI):

9681

UNIPROT ID:

O75140

Full Name:

DEP domain containing 5

Calculated MW:

1572 aa, 178 kDa

Purification Method:

Protein A purification

CloneNo.:

1F6G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51230-1, DEPDC5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60832-1-PBS. Detection antibody: 60832-2-PBS. Standard:Ag30176. Range: 0.781-100 ng/mL.