

For Research Use Only

MUM1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60834-1-PBS



Basic Information

Catalog Number:

60834-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29894

GenBank Accession Number:

BC008098

GeneID (NCBI):

84939

UNIPROT ID:

Q2TAK8

Full Name:

melanoma associated antigen (mutated) 1

Calculated MW:

18 kDa, 79 kDa

Purification Method:

Protein G purification

CloneNo.:

3E2C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51233-1, MUM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60834-1-PBS. Detection antibody: 60353-2-PBS. Standard:Ag3371. Range: 0.781-100 ng/mL.