

For Research Use Only

TTK Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60847-2-PBS



Basic Information

Catalog Number:

60847-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20556

GenBank Accession Number:

BC000633

GeneID (NCBI):

7272

UNIPROT ID:

P33981

Full Name:

TTK protein kinase

Calculated MW:

97 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E1E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51256-1, TTK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60847-1-PBS. Detection antibody: 60847-2-PBS. Standard:Ag20556. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51256-2, TTK Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60847-3-PBS. Detection antibody: 60847-2-PBS. Standard:Ag20556. Range: 0.195-100 ng/mL