

For Research Use Only

DUB3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60851-1-PBS



Basic Information

Catalog Number:

60851-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG24012

GenBank Accession Number:

BC157073

GeneID (NCBI):

377630

UNIPROT ID:

Q6R6M4

Full Name:

deubiquitinating enzyme 3

Purification Method:

Protein G Magarose purification

CloneNo.:

1G4E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51261-1, DUB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60851-1-PBS. Detection antibody: 60851-2-PBS. Standard:Ag24012. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP51261-2, DUB3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60851-1-PBS. Detection antibody: 60851-3-PBS. Standard:Ag24012. Range: 0.781-100 ng/mL