

For Research Use Only

GAS2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60858-2-PBS



Basic Information

Catalog Number:

60858-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29486

GenBank Accession Number:

BC040470

GeneID (NCBI):

2620

UNIPROT ID:

O43903

Full Name:

growth arrest-specific 2

Calculated MW:

313 aa, 35 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B3C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51271-1, GAS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60858-1-PBS. Detection antibody: 60858-2-PBS. Standard:Ag29486. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51271-4, GAS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60858-5-PBS. Detection antibody: 60858-2-PBS. Standard:Ag29486. Range: 0.781-100 ng/mL