

For Research Use Only

TNFRSF18 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60863-1-PBS



Basic Information

Catalog Number:

60863-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgM

Immunogen Catalog Number:

AG26103

GenBank Accession Number:

NM_004195

GeneID (NCBI):

8784

UNIPROT ID:

Q9Y5U5

Full Name:

tumor necrosis factor receptor
superfamily, member 18

Calculated MW:

26 kDa

Purification Method:

Thiophilic affinity chromatograph

CloneNo.:

1B1A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51275-1, GITR/TNFRSF18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60863-1-PBS. Detection antibody: 60863-2-PBS. Standard:Ag26103. Range: 6.25-100 ng/mL.