

For Research Use Only

FAM111B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60865-2-PBS



Basic Information

Catalog Number:

60865-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29380

GenBank Accession Number:

BC130539

GeneID (NCBI):

374393

UNIPROT ID:

Q6SJ93

Full Name:

family with sequence similarity 111, member B

Purification Method:

Protein G Magarose purification

CloneNo.:

1A8D11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51276-1, FAM111B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60865-1-PBS. Detection antibody: 60865-2-PBS. Standard:Ag29380. Range: 0.391-200 ng/mL