

For Research Use Only

IL22RA1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60866-1-PBS



Basic Information

Catalog Number:

60866-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29983

GenBank Accession Number:

BC029273

GeneID (NCBI):

58985

UNIPROT ID:

Q8N6P7

Full Name:

interleukin 22 receptor, alpha 1

Calculated MW:

574 aa, 63 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2A7B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51277-1, IL22RA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60866-1-PBS. Detection antibody: 60866-2-PBS. Standard:Ag29983. Range: 0.195-200 ng/mL



Cytometric bead array standard curve of MP51277-2, IL22RA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60866-1-PBS. Detection antibody: 60866-3-PBS. Standard:Ag29983. Range: 0.195-200 ng/mL