

For Research Use Only

# IL22RA1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 60866-2-PBS



## Basic Information

Catalog Number: 60866-2-PBS	GenBank Accession Number: BC029273	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 58985	CloneNo.: 1F8H3
Source: Mouse	UNIPROT ID: Q8N6P7	
Isotype: IgG1	Full Name: interleukin 22 receptor, alpha 1	
Immunogen Catalog Number: AG29983	Calculated MW: 574 aa, 63 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51277-1, IL22RA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60866-1-PBS. Detection antibody: 60866-2-PBS. Standard:Ag29983. Range: 0.195-200 ng/mL



Cytometric bead array standard curve of MP51277-3, IL22RA1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60866-2-PBS. Detection antibody: 60866-4-PBS. Standard:Ag29983. Range: 3.125-200 ng/mL