

For Research Use Only

H3F3B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60867-1-PBS



Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60867-1-PBS	BC001124	Protein G purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	3021	2B2A4
Source:	UNIPROT ID:	
Mouse	P84243	
Isotype:	Full Name:	
IgG1	H3 histone, family 3B (H3.3B)	
Immunogen Catalog Number:	Calculated MW:	
AG21191	15 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51278-1, H3F3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60867-1-PBS. Detection antibody: 60867-2-PBS. Standard:Ag21191. Range: 0.781-100 ng/mL.