

For Research Use Only

# H3F3B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60867-2-PBS



## Basic Information

Catalog Number: 60867-2-PBS	GenBank Accession Number: BC001124	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 3021	CloneNo.: 2A10B12
Source: Mouse	UNIPROT ID: P84243	
Isotype: IgG1	Full Name: H3 histone, family 3B (H3.3B)	
Immunogen Catalog Number: AG21191	Calculated MW: 15 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51278-1, H3F3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60867-1-PBS. Detection antibody: 60867-2-PBS. Standard:Ag21191. Range: 0.781-100 ng/mL.