

For Research Use Only

H3F3B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60867-4-PBS



Basic Information

Catalog Number:

60867-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21191

GenBank Accession Number:

BC001124

GeneID (NCBI):

3021

UNIPROT ID:

P84243

Full Name:

H3 histone, family 3B (H3.3B)

Calculated MW:

15 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D8A4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51278-2, H3F3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60867-3-PBS. Detection antibody: 60867-4-PBS. Standard:Ag21191. Range: 0.781-100 ng/mL.