

For Research Use Only

PCSK2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60870-1-PBS



Basic Information

Catalog Number:

60870-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG27919

GenBank Accession Number:

BC005815

GeneID (NCBI):

5126

UNIPROT ID:

P16519

Full Name:

proprotein convertase
subtilisin/kexin type 2

Calculated MW:

70.5 kDa

Purification Method:

Protein A purification

CloneNo.:

2G4C11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51284-1, PCSK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60870-1-PBS. Detection antibody: 60870-2-PBS. Standard:Ag27919. Range: 0.781-100 ng/mL.