

For Research Use Only

LGR5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60871-1-PBS



Basic Information

Catalog Number:

60871-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG16389

GenBank Accession Number:

BC096324

GeneID (NCBI):

8549

UNIPROT ID:

O75473

Full Name:

leucine-rich repeat-containing G protein-coupled receptor 5

Calculated MW:

907 aa, 100 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1G10B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51285-1, LGR5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60871-1-PBS. Detection antibody: 60871-2-PBS. Standard:Ag16389. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP51285-2, LGR5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60871-1-PBS. Detection antibody: 60871-3-PBS. Standard:Ag16389. Range: 0.195-100 ng/mL