

For Research Use Only

# DLX1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60872-2-PBS



## Basic Information

Catalog Number:	GenBank Accession Number:	Purification Method:
60872-2-PBS	BC036189	Protein G Magarose purification
Concentration:	GeneID (NCBI):	CloneNo.:
1 mg/ml	1745	1C10D9
Source:	UNIPROT ID:	
Mouse	P56177	
Isotype:	Full Name:	
IgG1	distal-less homeobox 1	
Immunogen Catalog Number:	Calculated MW:	
AG28661	255 aa, 27 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51286-1, DLX1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60872-1-PBS. Detection antibody: 60872-2-PBS. Standard:Ag28661. Range: 0.098-100 ng/mL.