

For Research Use Only

FTH1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60875-1-PBS



Basic Information

Catalog Number: 60875-1-PBS	GenBank Accession Number: BC000857	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 2495	CloneNo.: 4D3B2
Source: Mouse	UNIPROT ID: P02794	
Isotype: IgG1	Full Name: ferritin, heavy polypeptide 1	
Immunogen Catalog Number: AG21386	Calculated MW: 183 aa, 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51289-1, FTH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60875-1-PBS. Detection antibody: 60875-2-PBS. Standard:Ag21386. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51289-2, FTH1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60875-1-PBS. Detection antibody: 60875-3-PBS. Standard:Ag21386. Range: 0.781-50 ng/mL