

For Research Use Only

SLC10A1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60876-4-PBS



Basic Information

Catalog Number:

60876-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25680

GenBank Accession Number:

BC074724

GeneID (NCBI):

6554

UNIPROT ID:

Q14973

Full Name:

solute carrier family 10 (sodium/bile acid cotransporter family), member 1

Calculated MW:

349 aa, 38 kDa

Purification Method:

Protein G purification

CloneNo.:

1A11F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51290-3, SLC10A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60876-1-PBS. Detection antibody: 60876-4-PBS. Standard:Ag25680. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP51290-4, SLC10A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60876-5-PBS. Detection antibody: 60876-4-PBS. Standard:Ag25680. Range: 0.781-50 ng/mL.