

For Research Use Only

# CPEB1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60878-1-PBS



## Basic Information

Catalog Number:

60878-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28656

GenBank Accession Number:

BC035348

GeneID (NCBI):

64506

UNIPROT ID:

Q9BZB8

Full Name:

cytoplasmic polyadenylation  
element binding protein 1

Calculated MW:

486 aa, 54 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2E4E1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51293-1, CPEB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60878-1-PBS. Detection antibody: 60878-2-PBS. Standard:Ag28656. Range: 3.125-200 ng/mL.