

For Research Use Only

LILRB5 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60879-1-PBS



Basic Information

Catalog Number:

60879-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28624

GenBank Accession Number:

BC025704

GeneID (NCBI):

10990

UNIPROT ID:

O75023

Full Name:

leukocyte immunoglobulin-like receptor, subfamily B (with TM and ITIM domains), member 5

Calculated MW:

590 aa, 64 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G7D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51294-1, LILRB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60879-1-PBS. Detection antibody: 60879-2-PBS. Standard:Ag28624. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP51294-3, LILRB5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60879-1-PBS. Detection antibody: 60879-4-PBS. Standard:Ag28624. Range: 0.391-100 ng/mL