

For Research Use Only

# SPATS2L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60881-2-PBS



## Basic Information

Catalog Number:

60881-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9647

GenBank Accession Number:

BC018736

GeneID (NCBI):

26010

UNIPROT ID:

Q9NUQ6

Full Name:

viral DNA polymerase-transactivated protein 6

Calculated MW:

488 aa, 54 kDa

Purification Method:

Protein G purification

CloneNo.:

1E7A8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51296-1, SPATS2L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60881-1-PBS. Detection antibody: 60881-2-PBS. Standard:Ag9647. Range: 0.391-100 ng/mL.