

For Research Use Only

# PIK3C2B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60885-3-PBS



## Basic Information

Catalog Number:

60885-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20478

GenBank Accession Number:

BC132876

GeneID (NCBI):

5287

UNIPROT ID:

O00750

Full Name:

phosphoinositide-3-kinase, class 2, beta polypeptide

Calculated MW:

1634 aa, 185 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3E6F10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP51299-2, PIK3C2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60885-1-PBS. Detection antibody: 60885-3-PBS. Standard:Ag20478. Range: 0.781-100 ng/mL.