

For Research Use Only

PIK3C2B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60885-4-PBS



Basic Information

Catalog Number:

60885-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20478

GenBank Accession Number:

BC132876

GeneID (NCBI):

5287

UNIPROT ID:

O00750

Full Name:

phosphoinositide-3-kinase, class 2, beta polypeptide

Calculated MW:

1634 aa, 185 kDa

Purification Method:

Protein G purification

CloneNo.:

3D10C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51299-3, PIK3C2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60885-4-PBS. Detection antibody: 60885-5-PBS. Standard:Ag20478. Range: 0.781-100 ng/mL.