

For Research Use Only

GREM1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60888-3-PBS



Basic Information

Catalog Number:

60888-3-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12734

GenBank Accession Number:

BC093778

GeneID (NCBI):

26585

UNIPROT ID:

O60565

Full Name:

gremlin 1, cysteine knot superfamily, homolog (Xenopus laevis)

Calculated MW:

184 aa, 21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C1G2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51302-2, GREM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60888-1-PBS. Detection antibody: 60888-3-PBS. Standard:Ag12734. Range: 0.195-100 ng/mL.